

Calcium Alginate Dressing

2"x2" (5 x 5cm) 10/box GEN-13200	4"x4" (10 x 10cm) 10/box GEN-13500	5"x5" (12.5 x 12.5cm) 10/box GEN-13600	6"x6" (15.2 x 15.2cm) 10/box GEN-13700	10"x10" (25 x 25cm) 10/box GEN-13010
---------------------------------------------------------	-----------------------------------------------------------	---------------------------------------------------------------	---------------------------------------------------------------	-------------------------------------------------------------

- **Uniformly absorbs, collects and contains 20 times its weight in exudate**
- **Can be molded and folded to fit the size of the wound**
- **Easy to remove – does not adhere to the healing tissue of the wound**



Gentell's Calcium Alginate Dressings are a sterile, comfortable, advanced fiber-structured alginate with a highly absorbent capacity. Alginate dressings absorb, collect, and contain exudate while providing a moist healing environment. A reaction between the calcium in the dressing and the sodium in the wound exudate creates a gel-like substance that promotes moist wound healing. Use on leg ulcers, pressure ulcers, diabetic foot ulcers, and surgical wounds. May be used for minor conditions such as lacerations, abrasions, skin tears, and minor burns.

Directions

1. Irrigate the wound with Gentell Dermal Wound Cleanser and gently dry the skin surrounding the wound site.
2. Choose the appropriate size dressing based on the dimensions of the wound.
3. Cover wound with a secondary dressing such as Gentell Super Absorbent, Gentell Comfortell™ or Gentell Bordered Foam Dressing.
4. Repeat daily or as ordered by a physician.